

this  
**Webinar** is powered by  
**PerkinElmer**

**23 June 2022**

3:00 pm – 4:00 pm | CEST, Berlin  
5:00 pm – 6:00 pm | GST, Dubai  
9:00 am – 10:00 am | EDT, New York City  
6:00 am – 7:00 am | PDT, Los Angeles

pV magazine  
**webinars**

# Chasing the sun – advanced analysis for innovative PV materials



**Tim Sylvia**

Editor  
pv magazine USA



**Lucy Jenner**

Segment Manager of Industrial  
**PerkinElmer**



**Dr. Gernot Oreski**

Division Manager for  
Sustainable Polymer Solutions  
**Polymer Competence Center  
Leoben**



**Peter van Nijnatten**

R&D Director  
**OMT Solutions**



**Krystelle Mafina**

UVS/TEA Product Manager  
**PerkinElmer**

# Welcome!

Do you have any questions? ? 

Send them in via the Q&A tab.  We aim to answer as many as we can today!

You can also let us know of any tech problems there.

We are recording this webinar today. 

We'll let you know by email where to find it and the slide deck, so you can re-watch it at your convenience.  

# Chasing the Sun – Advanced Analysis for Innovative PV Materials


Advanced Analysis for Innovative PV Materials

Gerlinde Wita, Gernot Oreski, Serge Timmermans, Krystelle Mafina




# Solar Economics & The Future

Costs continue to fall for solar and wind power technologies



- Fastest growing of the renewables
- Huge diversity of manufacturers – facing similar challenges
- Cheapest since 2020 reporting





# Innovation & Material Science Challenges



- Material Innovation – esp glass, silicon, polymers
- Properties of cells, wafers, interconnections is vital to understand
- High product quality and reduce costs are key in order to stay ahead of the competition

# Polymers for photovoltaics


## Ensuring module reliability via proper material testing

**Dr. Gernot Oreski**

Polymer Competence Center Leoben GmbH,  
Roseggerstraße 12, 8700 Leoben  
+43 664 88679331  
gernot.oreski@pccl.at



- Continuous strong market growth [1]
  - ✓ *Annual growth rates between 30 and 40 %*
- PV market is highly competitive with enormous cost pressure
  - ✓ *Average learning rate of about 80%*
- Market is driven by technological developments and market conditions [1]
  - ✓ *In recent years many new PV module materials and PV module designs were launched on the market*
- Motivation for new PV module materials and module designs [2]
  - ✓ *Decrease of LCOE: Cost reduction and performance improvement*
  - ✓ *New technological requirements*
  - ✓ *Sustainability and legal regulations*




**Average price of PV modules fell by 20% for each doubling of production volume**

[1] European Commission, PV Status Report 2019, <https://ec.europa.eu/jrc/en/publication/eur-scientific-and-technical-research-reports/pv-status-report-2019>


[2] G Oreski et al 2022 Prog. Energy 4 032003, <https://www.doi.org/10.1088/2516-1083/ac6f3f>

Wafer Size Comparison M0 - M12+ © RENMA Technologies GmbH




<https://www.rena.com/en/products/large-wafer-wet-processing/>

## New materials and technologies in PV modules [2]



<https://doi.org/10.1016/j.apsusc.2020.145420>




Structured foil

<https://www.sunportpower.com/>

- Wafer size increases
- PERC and HJT will become dominant cell technologies
- Bifacial DG PV modules
- Polyolefin encapsulants
- New backsheet architectures (co-extruded, coated)
- New interconnection technologies
- Half cells will gain more market share
- Module sizes increase

### New IEA Task 13 report on DESIGNING NEW MATERIALS FOR PHOTOVOLTAICS

- ✓ Report IEA-PVPS T13-13:2021, April 2021
- ✓ <https://iea-pvps.org/key-topics/designing-new-materials-for-photovoltaics/>



<https://www.cleanenergyreviews.info/blog/most-efficient-solar-panels>



## PV systems designed for specific environmental conditions / Integrated PV

Desert PV



Vehicle integrated PV



© Sono Motors

Agri PV



© Fraunhofer ISE

BIPV



© Ertex Solar

PV in transport infrastructure



© AIT


Floating PV



© Baywa r.e.



- Check of compatibility of PV module components will get more and more important in the future, as the variety on materials and components will grow
- **Emergence of new degradation modes** (e.g. LeTid, PID.... )




## Adhesion - delamination

- ✓ Adhesion to glass and solar cell strongly dependent of lamination parameters
- ✓ Surface treatment of backsheets usually optimized for adhesion to EVA but not alternative encapsulants



**Backsheet yellowing:** Migration of additives into backsheet-encapsulant interface are main cause for backsheet yellowing




**Corrosion:** Broad variety of new ribbon materials, interconnection technologies and encapsulant films

**Completely new degradation mechanisms?**

**Constant need for adaption of test methods and standards**

# Why is it important to know your materials?

## *Example: Encapsulant degradation*




### **(6) EFFECTS OF ARTIFICIAL AGEING TESTS ON EVA DEGRADATION: INFLUENCE OF MICROCLIMATE AND METHODOLOGY APPROACH**

**Chiara Barretta<sup>1</sup>, Gernot Oreski<sup>1</sup>, Nikoleta Kyranaki<sup>2,3</sup>, Djamel E. Mansour<sup>4</sup>, Thomas R. Betts<sup>2</sup>, Luciana Pitta Bauermann<sup>4</sup>, Katharina Resch-Fauster<sup>5</sup>.**

<sup>1</sup>Polymer Competence Center Leoben GmbH (PCCL), Leoben, Austria <sup>2</sup>Centre for Renewable Energy Systems (CREST), Loughborough University, Loughborough, United Kingdom <sup>3</sup>Commissariat à l'Energie Atomique, CEA Grenoble, DRT/LITEN/DTS/SCPV/LCT, Le Bourget-du-Lac, France <sup>4</sup>Fraunhofer Institute for Solar Energy Systems (ISE), Freiburg im Breisgau, Germany <sup>5</sup>Department of Material Science and Testing of Plastics, University of Leoben, Leoben, Austria



## Additive analysis of EVA from PV modules




*Thermal Desorption - Gas chromatography  
coupled with mass spectrometry*

### Ageing tests:


- DH (85°C/85% R.H.), 3000 hours
- UV exposure ~85 kWh/m<sup>2</sup> (ISO 4892-3 Cycle 1)
- Single layers of cured encapsulant

- 1) Acetic acid
- 2) Fragment of peroxide
- 3) Antioxidant
- 4) Crosslinking agent
- 5) Adhesion promotor
- 6) UV absorber






Artificial ageing test	Main parameters	Maximum Time / Dose
<b>Damp Heat</b> (IEC 61215- 2:2016 MQT 13)	85 °C, 85% Relative Humidity (RH)	2000 h
<b>Dry UV</b> IEC 61215-2:2016 MQT 11	250 W m <sup>-2</sup> (280 nm – 400 nm), 60 °C	500 kWh m <sup>-2</sup>
<b>UV-DH combined</b>	180 W m <sup>-2</sup> , 60 °C, 85% RH	250 kWh m <sup>-2</sup>



- ✖ EVA above the cell
- ✖ EVA above the backsheet
- ✖ EVA outside

Different microclimate at each measurements position

- Decrease of  $P_{MAX}$  of about 5% for the sample exposed to UVDH combined test.
- Decrease of power output associated to decrease in  $I_{sc}$  with the same trend
- Decrease in  $I_{sc}$  might be associated to **discoloration** of encapsulant and/or **cell cracks**




<b>EVA above the backsheet</b>	<b>REF</b>	<b>DH</b>	<b>UV</b>	<b>UV-DH</b>
Antioxidant - BHT	✓	✓	n.d.	n.d.
UV absorber - Benzophenone	✓	✓	✓	✓
<b>EVA above the cell</b>	<b>REF</b>	<b>DH</b>	<b>UV</b>	<b>UV-DH</b>
Antioxidant - BHT	✓	✓	✓	n.d.
UV absorber - Benzophenone	✓	✓	✓	✓
<b>EVA outside</b>	<b>REF</b>	<b>DH</b>	<b>UV</b>	<b>UV-DH</b>
Antioxidant - BHT	✓	✓	n.d.	n.d.
UV absorber - Benzophenone	✓	✓	n.d.	✓

- The antioxidant is no longer detectable in the samples exposed to UV and UV-DH combined test
- The UV absorber is detectable for all samples

✓ = detected, n.d. = not detected


*Thermal Desorption - Gas chromatography coupled with mass spectrometry*



*UV/VIS/NIR spectroscopy*

- No relevant discoloration was observed for the EVA above the cell
- Discoloration was observed for the EVA above the backsheet, regardless the artificial ageing test applied
- Different mechanism behind discoloration
  - ✓ *DH: Interaction at backsheet / encapsulant interface*
  - ✓ *UV exposed samples: Discoloration correlated to degradation products of the antioxidant*


## Oxidation index



## FTIR-ATR spectroscopy

- No relevant changes detected in the encapsulant above the backsheet and above the cell, regardless the artificial test applied
- Strong oxidation signs for the encapsulant directly exposed to the environment in the UV and UV-DH combined test

## FTIR ATR spectroscopy of EVA outside




- Strong influence of micro-climate - no EVA degradation within the PV module
- Power loss can be attributed to yellowing
- Yellowing is not caused by degradation of the EVA molecules, but by additives

**How to predict and possibly prevent  
polymer related PV module failure modes?**

**Proper material selection  
and testing, testing, testing!**




## Quality assurance to increase lifetime and reliability



### Sequential / combined stress testing

- Currently the industry relies mostly on extended IEC testing for qualification of new module materials and module designs
- Material interactions and incompatibilities are getting in the focus of material and module developers
- Simultaneous combined or sequential stresses (UV, humidity, temperature and thermo-mechanical load) lead to more realistic degradation of PV modules in lab testing



### Testing of encapsulants & backsheets in the lab

- Different levels with increasing effort in sample preparation and test procedures
- All levels of testing needed in order to guarantee high quality and reliability

- **Role of the polymers in photovoltaic energy generation has generally been underestimated**
  - ✓ *No active role in power generation*
- **Choice of polymers has distinctive impact on PV modules attributes such as**
  - 1) ***Efficiency**, as the optical properties of encapsulant (transmittance) and backsheet (reflectance, back scattering) define the number of photons arriving at the solar cell*
  - 2) ***Quality**, as the main infant failures are caused by bad processing parameters, which are defined by the encapsulant properties, and material incompatibilities*
  - 3) ***Reliability**, as most PV module degradation modes are directly linked to polymer degradation and material interactions with polymer components*

→ **Better understanding of material properties of polymers in PV modules and their influence degradation processes is a precondition for a successful development of new components and reliable PV module designs**

## Thank you for your attention!

### Project funding



This project has received funding from the European Union's Horizon 2020 programme under GA. No. 721452.



Energy Research Programm -, FFG No. 867267, Klima- und Energiefonds



Energy Research Programm -, FFG No. 850414, Klima- und Energiefonds



<https://orcid.org/0000-0003-4223-9047>



# Optical Characterisation of PV Materials using the Lambda 1050

Peter van Nijnatten



# High-End UV(Vis/NIR) applications for PV R&D

- **Variable Angle Spectroscopy**
  - Thin Film Analysis (layer thickness, Optical Constants)
  - Angular Resolved Scattering for PV Light Management
- **Measuring PV Cover Glass**
  - <sup>n</sup> Glass Redox Measurements ( $\text{Fe}^{2+}/\text{Fe}^{3+}$ )
  - Total Transmittance and Reflectance measurements (diffuse and patterned cover glass)

<sup>k</sup>



# Tools for Variable Angle Spectroscopy

## ARTA






## TAMS





# Thin Film Analysis by VAS

Non-linear regression is used to minimise the difference between measured and calculated spectra by varying the values for the thickness of the layers and the parameters of the dispersion models, resulting in their best estimates.




# Angular Resolved Scattering





# Surface Topology of TCO's analysed by Angular Resolved Scattering



**Figure 2.** AID of light that is scattered when transmitting a surface-textured layer of  $\text{SnO}:\text{F}_2$  (a),  $\text{ZnO}:\text{Al}$  (b),  $\text{ZnO}:\text{B}$  (c) and a 1D grating (d). The AID is shown as a function of the scattering angle  $\theta$



# Surface Topology of TCO's analysed by Angular Resolved Scattering









Figure 9: Integrated angular resolved light scattering (left) for a very rough and steep morphology (right).

# Characterising Solar Cel cover glass





# 2007 Inter-Laboratory comparison of transmittance measurements







## Measuring Solar Cell cover glass with the standard 150 mm sphere



# UL270 Sphere accessory with fixed 8° AOI for Transmittance and Reflectance





# 2017 NFRC Intercomparison



**National  
Fenestration  
Rating Council  
Incorporated**


NFRC 300-2017[E0A2]

Test Method for  
Determining the Solar Optical Properties  
of Glazing Materials and Systems


© 2013 NATIONAL FENESTRATION RATING COUNCIL, INC.

**PREPARED BY:**

National Fenestration Rating Council  
6305 Ivy Lane, Suite 140  
Greenbelt, MD 20770-6323  
Voice: (301) 589-1776  
Fax: (301) 589-3884  
Email: [info@nfr.org](mailto:info@nfr.org)  
Website: [www.nfr.org](http://www.nfr.org)





# Measuring Directional R&T of diffuse and patterned cover glass





- Sample size up to 400 x 300 mm
- 270 mm integrating sphere with a large (100 mm) sample port
- PMT for UV/Vis, InGaAs detector for NIR
- Wavelength range 250 nm – 2,400 nm
- Measurement spot of approximately 10 x 25 mm at 0 degrees incidence
- Separate P and S polarization (automated polarizer drive)
- Directional Reflectance at 8° – 70° AOI (up to 80° with lower accuracy)
- Directional Transmittance at 0° – 70° AOI (up to 80° with lower accuracy)
- Haze measurement at -8° AOI in Transmittance

## Measuring Directional R&T of diffuse and patterned cover glass









## Measuring Directional R&T of diffuse and patterned cover glass







# Colored Cover Glass for Building Integrated PV



# Examples of Directional R&T Results








# Analytical Technology Pros & Cons

Krystelle Mafina

# Our Materials Characterisation Product Portfolio





# InfraRed for the Chemical Characterisation

Infrared spectroscopy (IR) measures vibrational and rotational transitions

Allows for both qualitative and quantitative information about functional groups

The range is split in 3 regions (FIR, MIR & NIR) each giving rise to a different set of information depending on need

Mostly used for chemical characterisation as it gives rise to accurate data when sampling is done appropriately



# Pros & Cons of Infrared Spectroscopy

BENEFITS	LIMITATIONS
MIR/NIR/FIR in one instrument available	Tri-range in one system can cause some loss in sensitivity
Meeting a wide range of samples and method development flexibility	Limited use of all accessories for some of the regions especially FIR
Enhanced ATR accessory available	Only measures surface layer and a small area, a few microns depth at most
Faster scanning rate for maximum scan per seconds	
Higher accuracy of library searching using ATR	Spectra different from the ones obtained via transmission, must create an extensive one
One vendor Hyphenation with single instrument EGA solution	Amount of sample to use and destructive, always begin the run with thermal analysis

# Hyphenation for the Chemical Characterisation

When at least 2 systems are combined to generate 2 sets of data from 1 sampling preparation

Usually based on the analysis of gaseous analytes for a better sampling transfer between the systems

Heat is involved thus care must be taken into maintaining heated sample transfer from one system to another

Carrier gas is required to enable a smooth sampling transfer

Event triggers are needed to start or stop the various system collecting data at the time/weight/temperature required for the analysis



# Pros & Cons of Hyphenation




Instrument	TG-IR	TG-GCMS	TG-IR-GCMS
Consumables	Pans (crucibles) carrier gas (He/N <sub>2</sub> )	Pans (crucibles), columns, carrier gas(He)	Pans (crucibles), columns, carrier gas (He)
Software	Spectrum IR	Spectrum IR and TurboMass	Spectrum IR and Turbomass
BENEFITS	<ul style="list-style-type: none"> <li>Resulting data collected are time base resolved for a qualitative interpretation</li> <li>Ready in seconds - no alignment required</li> </ul>	<ul style="list-style-type: none"> <li>Unique injection valve system Allow the user to do TGA-GCMS or TGA-MS</li> <li>Snap-shot or real-time measurement</li> </ul>	<ul style="list-style-type: none"> <li>Highest controlled temperature transfer line to prevent condensation (Analysis of high boiling point samples)</li> </ul>
LIMITATIONS	Quantification is available but not straight-forward	Cross-contamination within column, thus maintenance is important	Cross-contamination within column, thus maintenance is important





# Our Materials Characterisation Product Portfolio


**Chemical Characterization**


Raw Materials Identification  
Polymer/Material Characterisation  
Product defects & inclusions

**PolymerID Analyzer**


**FTIR/FTNIR**


**Microscopy/Imaging**


**TG-IR**

**TG-IR-GC/MS**

Residual monomer/polymer content  
Extractables & Leachables, Fogging, PCB/PAH  
VOC, Phtalate Analysis, Softener, Bisphenol A


**GC & GC/MS**


**Headspace**


**LC & LC/MS**


**Mechanical & Thermal Characterisation**


Mechanical & Thermal Characterisation


**TG-IR-GC-MS**

**DSC**


**TGA**


**STA**


**DMA**


**TMA**


Analytics of Additives  
Elemental Analysis  
Colour Measurements


**ICP-OES**


**LC & LC/MS**


**UV-Vis**

**High-End UV-Vis**

**OneSource**  
Laboratory Services

**TIBCO Spotfire**

**NetPlus**

**E-Notebook**

41

41

# UV-Vis-NIR for Mechanical Characterisation


Exists in 2 principles

1. PDA design
2. Double Beam design

Photodiode array spectrophotometers work differently from conventional spectrophotometers as the sample is irradiated with white light and then the light passing through the sample is dispersed (using a polychromator) onto a photodiode detector.

Detection is simultaneous at all wavelengths and consequently the time taken to acquire a complete full range spectrum is less than two seconds.

Conventional instruments, by comparison, irradiate the sample one wavelength at a time (i.e. they are sequential).



# Pros & Cons of UV-Vis-IR

BENEFITS	LIMITATIONS
Double beam instruments collect data sequentially for better precision and accuracy	The data collection is slow as it does one wavelength at one time
Double beam instrument can handle higher absorbance	
Double instruments are complex system with moving parts	Require good maintenance
Covers a full UV-Vis-NIR range	PDA only covers UV-Vis
PDA instruments have fast scanning	Loss of sensitivity
PDA system have no moving parts and no need for sample lid	
Sample compartment on the double beam instrument is large for a variety of sample size	May require a specific accessory

# Thermal Analysis Product Portfolio

## Thermal Characterization in QA/QC

Raw Materials Identification  
Polymer/Material Characterisation



Single furnace DSC  
(Heatflux)



TGA for routine  
applications



Simultaneous TG/DTA  
(DSC)



PerkinElmer<sup>®</sup>  
*For the Better*



Software Pyris<sup>™</sup>

## Mechanical & Thermal Characterisation for Research & problem solving

Thermal Characterisation



Double furnace DSC  
(Power Compensation)

TGA for High Demand  
& Hyphenation

Thermal Mechanical Characterisation



TMA & DMA for various applications

HYPHENATION - Get better insights to understand when something happened, and also what happened





# Thermal Analysis for Q&A

These systems are mostly used for routine analysis, especially useful for characterising the polymeric materials (layer) in the solar cells


They are used to identify and help discriminate between polymeric materials and their resulting effects

They are robust

Allow for high temperature analysis up to 1600 deg

They measure physical properties of a material against temperature or time.

- Glass transition
- Decomposition
- Degree of crystallisation
- Melting point etc ...



DSC6000/4000



TGA4000



STA8000/6000

# Thermal Analysis for Problem Solving & Research

These instruments are more sensitive, allowing for more in-depth understanding

More accurate balance in the TGA8000

More parameters to be collated in the DMA8000

During these analyses, sample preparation is key for optimum results

They measure mechanical properties of a material against temperature or time.

- Mechanical behaviour
- Weight loss
- Long-term behaviour
- Cross-linking properties etc ...



DMA8000

# Pros & Cons of Thermal Analysis

BENEFITS	LIMITATIONS
A variety of properties can be measures using thermal product line	Care must be taken when choosing the system for analysis
Up to 1.2 g of sample can be analysed	Only small density of sample can be used
Small density of material equates to improved precision and accuracy	May not be representative
Identify the material's minimum process temperature	May not be able to differentiate if the material is too similar
Identify the amount of energy required to melt the material (heat of fusion)	
Different carrier gas can be used to mimic environmental conditions	Only up to 3 gases can be used .. but there is mixing pod accessory available
Identify, quantify and characterise material's components	Often struggles with multi-component dependent on sample preparation



# Renewable Energy & Sustainability

## Our Commitment to Sustainability

Our dedication to a more sustainable world can be seen in our three pillars: instruments and products, customer and industry focus, and internal operations.

By incorporating thoughtful action in each of these areas, we enable businesses like yours to successfully achieve sustainability objectives while minimizing the environmental impact.

Instruments  
and Solutions

Sustainability by  
Design


Industry and  
Customer Focus

Enabling  
Sustainable  
Innovation

Internal  
Operations  
Sustainable  
Goals and  
Values

EMPOWERING  
SUSTAINABILITY™

[www.perkinelmer.com/category/sustainability-solutions](http://www.perkinelmer.com/category/sustainability-solutions)



Thank You.  
Any Questions?



For More Information  
[perkinelmer.com/category/solar](https://perkinelmer.com/category/solar)



this  
**Webinar** is powered by  
**PerkinElmer**

**23 June 2022**

3:00 pm – 4:00 pm | CEST, Berlin  
5:00 pm – 6:00 pm | GST, Dubai  
9:00 am – 10:00 am | EDT, New York City  
6:00 am – 7:00 am | PDT, Los Angeles

pV magazine  
**webinars**

## Chasing the sun – advanced analysis for innovative PV materials

### Q&A



**Tim Sylvia**

Editor  
pv magazine USA



**Lucy Jenner**

Segment Manager of Industrial  
**PerkinElmer**



**Dr. Gernot Oreski**

Division Manager for  
Sustainable Polymer Solutions  
**Polymer Competence Center  
Leoben**



**Peter van Nijnatten**

R&D Director  
**OMT Solutions**



**Krystelle Mafina**

UVS/TEA Product Manager  
**PerkinElmer**



# The latest news | print & online



**10% off**  
your subscription  
with  
**Webinars10**



## Half a penny for “near firm” solar and trillions in renewable opportunities

by John Fitzgerald Weaver



Most-  
read  
online!

## Biden halts solar tariffs for two years. What’s next?

by Ryan Kennedy and Tim Sylvia



# Coming up next...

## Thursday, 30 June 2022

5:00 pm – 6:00 pm CEST, Berlin, Madrid

11:00 am – 12:00 pm EDT, New York City

## Tuesday, 5 July 2022

11:00 am – 12:00 pm CEST, Berlin, Madrid

10:00 am – 11:00 am GMT, London

**Many more to come!**

**Alleviating  
battery storage  
growing pains  
through quality  
assurance**

**The pressing need for  
powerful storage:  
Technology trends,  
safety optimization,  
and challenges for  
deployment**

In the next weeks, we will continuously add further webinars with innovative partners and the latest topics.

Check out our pv magazine Webinar program at:

[www.pv-magazine.com/webinars](http://www.pv-magazine.com/webinars)

Registration, downloads  
& recordings are also be  
found there.



this  
**webinar** is powered by  
PerkinElmer

pv magazine  
**webinars**



**Tim Sylvia**  
Editor  
pv magazine USA

**Thank you for  
joining today!**